


**BE1431 Visit telephone transmitter**

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**36** User manual

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**68** Käyttöopas

**84** Mode d'emploi

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**116** Brukermanual

**132** Bruksanvisning

**149** Appendix (EN)

DE

DK

EN

ES

FI

FR

NL

NO

SE

App

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
DK

## Read this first

This booklet guides you on how to use and maintain your new transmitter. Ensure you read this booklet carefully, including the **Warnings** section. This will help you get the most out of your product. If you have additional questions, contact your hearing care professional.


|                               |                                                                                                                                                                                                                                                                                                                                                |
|-------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Intended purpose</b>       | The intended purpose of the Visit system is to alert deaf and hard of hearing people of important signals in their home.                                                                                                                                                                                                                       |
| <b>Intended user</b>          | A person with mild to severe hearing loss or deafness.                                                                                                                                                                                                                                                                                         |
| <b>Intended user group</b>    | The intended user group consists of people of all ages with mild to severe hearing loss or deafness that need audio, visual or sensory amplification.                                                                                                                                                                                          |
| <b>Principle of operation</b> | The Visit alerting system consists of a set of wirelessly connected transmitters and receivers that are located across the home. When a transmitter detects an activity, it signals the Visit receiver that alerts the user with sound, flashes or vibrations. In order to use the system, you need at least one transmitter and one receiver. |

## Overview


## Installation

- 1 Remove the battery pull tab to start the unit. Clean the wall with the wet wipe and remove the protective film from the Velcro.
- 2 Mount the transmitter on the wall. You can also use the supplied screw and plug.
- 3 Connect the telephone splitter between the telephone plugs and the socket.

## Testing

- 1 Use your mobile phone or ask a friend to call your landline telephone.
- 2 The yellow Visit LED on the receiver lights up and it starts to sound, flash or vibrate, depending on the receiver. If it doesn't react, read **Troubleshooting**.

## General warnings

This section contains important information about safety, handling and operating conditions. Keep this booklet for future use. If you are just installing the device, this booklet must be given to the householder.

### Hazard warnings

- Failure to follow these safety instructions could result in fire, electric shock, or other injury or damage to the device or other property.
- Make sure to keep this device out of reach of children under 3 years.
- Note that alarms and notifications can be missed if batteries run out.
- Use only battery types that are specified in this booklet.
- This device is designed for indoor use only. Do not expose the device to moisture.
- Protect the device from shocks during storage and transport.
- This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the user guide, may cause harmful interference to radio communication.

- Do not use or store this device near any heat sources such as naked flames, radiators, ovens or other devices that produce heat.
- Do not dismantle the device; there is a risk of electric shock. Tampering with or dismantling the device will void warranty.
- Protect cables from any potential source of damage.
- Do not make any changes or modifications to this device. Use only original Bellman & Symfon accessories to avoid any electrical shock.
- Batteries are toxic. Do not swallow them! Keep out of reach of children and pets. If they are swallowed, consult your physician immediately!

### Information on product safety

- Failure to follow these instructions could result in damage to the device and void the warranty.
- Do not expose batteries to fire or to direct sunlight.
- Do not drop your device. Dropping onto a hard surface can damage it
- Do not use the device in areas where electronic equipment is prohibited.

- The device may only be repaired by an authorized service center.
- If a serious incident occurs in relation to this device, contact the manufacturer and relevant authority.
- If you encounter other problems with your device, contact the point of purchase, your local Bellman & Symfon office or the manufacturer. Visit **bellman.com** for contact info.
- Without permission granted by the NCC, no company, enterprise or user is allowed to change the frequency, enhance the transmitting power or alter the original characteristics or performance on an approved low power radio frequency device.
- This device should not be used on aircrafts unless specifically permitted by flight personnel.

## Operating conditions

Operate the device in a dry environment within the temperature and humidity limits stated in this booklet. If the product gets wet or is exposed to moisture, it should no longer be regarded as reliable and should therefore be replaced.

## Cleaning

Disconnect all cables before you clean your device. Use a soft, lint-free cloth. Avoid getting moisture in openings. Do not use household cleaners, aerosol sprays, solvents, alcohol, ammonia, abrasives or compressed air.

## Service and support

If the product appears to be damaged or doesn't function properly, follow the instructions in this booklet. If the product still doesn't function as intended, contact your local hearing care professional for information on service and warranty.

## Warranty conditions

Bellman & Symfon guarantees this product (excluding the battery) for two (2) years from date of purchase against any defects that are due to faulty materials or workmanship. This guarantee only applies to normal conditions of use and service, and does not include damage resulting from accident, neglect, misuse, unauthorized dismantling, or contamination howsoever caused. This guarantee excludes incidental and consequential damage. Further the warranty does not cover Acts of God, such as

fire, flood, hurricanes and tornadoes. This warranty gives you specific legal rights and you may also have other rights that vary from state to state. Some states or jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitations on how long an implied warranty last so the above limitation may not apply to you. This guarantee is in addition to your statutory rights as a consumer. The above warranty may not be altered except in writing signed by both parties hereto.

### Configuration options

This transmitter can be configured with the following Visit receivers and accessories:

**Compatible Visit receivers**

- BE1450 Visit portable receiver
- BE1441 Visit flash receiver
- BE1442 Visit flash receiver
- BE1580 Visit alarm clock receiver
- BE1470 Visit pager receiver
- BE1560 Visit wrist receiver

**Available accessories**

- BE9026 Contact mat
- BE9023 Magnetic switch
- BE9253 External trig cable

For detailed product information, see the corresponding user manual.

### Compliance information

Hereby Bellman & Symfon declares that, in Europe, this product is in compliance with the essential requirements of the Medical Device Regulation EU 2017/745 as well as the directives and regulations listed below. The full text of the declaration of conformity can be obtained from Bellman & Symfon or your local Bellman & Symfon representative. Visit **bellman.com** for contact information.

Radio Equipment Directive (RED)

Medical Device Regulation (MDR)

Waste Electrical & Electronic Equipment (WEEE)

RoHS Directive

REACH Regulation

EC Battery Directive

## Regulatory symbols

|                                                                                    |                                                                                                                                                         |
|------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | With this symbol, Bellman & Symfon confirms that the product meets the Medical Device Regulation EU 2017/745.                                           |
|    | This symbol indicates the manufacturer's serial number so that a specific medical device can be identified. It's available on the product and gift box. |
|    | This symbol indicates the manufacturer's catalogue number so that the medical device can be identified. It's available on the product and gift box.     |
|    | This symbol indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.                                  |
|  | This symbol indicates that the user should consult the instruction for use.                                                                             |

|                                                                                     |                                                                                                                                                |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
|  | This symbol indicates that it is important for the user to pay attention to the relevant warning notices in the user guide.                    |
|  | This symbol indicates important information for handling and product safety.                                                                   |
|  | Temperature during transport and storage: -10° to 50° C, 14° – 122° F<br>Temperature during operation: 0° to 35° C, 32° to 95° F               |
|  | Humidity during transportation and storage: <90%, non-condensing<br>Humidity during operation: 15% – 90%, non-condensing                       |
|  | Atmospheric pressure during operation, transportation and storage:<br>700 hpa to 1060 hpa                                                      |
| <b>Operating conditions</b>                                                         | This device is designed such that it functions without problems or restrictions if used as intended, unless otherwise noted in the user guide. |

With this CE symbol, Bellman & Symfon confirms that the product meets EU standards for health, safety, and environmental protection as well as the Radio Equipment Directive 2014/53/EU.

This symbol indicates that the product shall not be treated as household waste. Please hand over your old or unused product to the applicable collection point for the recycling of electrical and electronic equipment or bring your old product to your hearing care professional for appropriate disposal.

## ISO Certification of legal manufacturer

Bellman is certified in accordance with SS-EN ISO 9001 and SS-EN ISO 13485.  
SS-EN ISO 9001 Cert. No.: CN19/42071, SS-EN ISO 13485 Cert. No.: CN19/42070

### Certification Body

SGS United Kingdom Ltd, Rossmore Business Park Ellesmere Port Cheshire CH65 3EN UK

## Technical specifications

|                          |                                                                                                                                                                         |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Dimensions</b>        | 100 × 65 × 27 mm, 4.0 × 2.6 × 1.1"                                                                                                                                      |
| <b>Weight</b>            | 134 g, 4.7 oz. including battery and telephone cord                                                                                                                     |
| <b>Power and battery</b> | Batteries: 2 × 1.5 V AA (LR6) lithium or alkaline<br>Power consumption: Active < 70 mA, Idle position < 15 µA<br>Operating time: Alkaline ~ 5 years, Lithium ~ 10 years |
| <b>Activation</b>        | Via test buttons, RJ11 input to analogue telephone network (26-120V RMS, 15-100Hz), 2.5 mm ext trig, 3.5 mm ext trig                                                    |
| <b>Max. output power</b> | 4.7 dBm                                                                                                                                                                 |
| <b>Radio frequency</b>   | 868.30MHz                                                                                                                                                               |
| <b>Coverage</b>          | 50 - 250 m depending on the characteristics of the building.                                                                                                            |
| <b>No. of radio keys</b> | 64 as standard. Up to 256 with optional software.                                                                                                                       |
| <b>In the box</b>        | <ul style="list-style-type: none"><li>▪ BE1431 Visit telephone transmitter</li><li>▪ 1.5V AA alkaline batteries</li><li>▪ Velcro, screws and plugs</li></ul>            |

EN

## Troubleshooting

| If                                                                                          | Try this                                                                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The LEDs blink in orange every minute                                                       | <ul style="list-style-type: none"><li>▪ The batteries are depleted. Replace them with 1.5 V AA (LR6) lithium or alkaline batteries.</li></ul>                                                                                                                                                                                              |
| The transmitter LED lights up in green but the receiver doesn't respond                     | <ul style="list-style-type: none"><li>▪ Check the transmitter batteries and the receiver batteries and connections.</li><li>▪ Move the receiver closer to the transmitter to make sure it's within radio range.</li><li>▪ Check that the units are set to the same radio key, see <b>Changing the radio key</b> in the Appendix.</li></ul> |
| The transmitter LED doesn't light up when the phone rings or when an accessory is triggered | <ul style="list-style-type: none"><li>▪ Press the test buttons on the transmitter. If the LED lights up in green, check all connections. If the LED doesn't light up in green, replace the batteries. Only use 1.5 V AA (LR6) lithium or alkaline batteries.</li></ul>                                                                     |

| If                                                                 | Try this                                                                                                                                                                                                                   |
|--------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The transmitter LED doesn't light up when I press the test buttons | <ul style="list-style-type: none"><li>▪ Replace the batteries. Only use 1.5 V AA (LR6) lithium or alkaline batteries. If the LED still doesn't light up, contact your retailer for service information.</li></ul>          |
| The receiver is activated for no apparent reason                   | <ul style="list-style-type: none"><li>▪ There is probably another Visit system installed nearby that triggers your system. Change the radio key on all units, see <b>Changing the radio key</b> in the Appendix.</li></ul> |

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
## Testing the signal

### Using the test buttons

- 1 Press both test buttons simultaneously on the transmitter. The top LED lights up in green to show that a radio signal is being transmitted.
- 2 The yellow Visit LED on the receiver lights up to show that the signal was received. In addition, it starts to sound, flash or vibrate depending on the receiver.

**Note:** If you have changed the signal pattern, it will react in accordance with the table on the following spread.


## Testing the signal

### Using a mobile phone

- 1 Use for instance a mobile phone to call the landline telephone. The top LED on the transmitter lights up in green to show that an incoming call is detected.
- 2 The yellow Visit LED on the receiver lights up to show that the radio signal was received. In addition, it starts to sound, flash or vibrate depending on the receiver.

**Note:** If you have changed the signal pattern, it will react in accordance with the table on the following spread.

## Default signal pattern

When the telephone transmitter is activated by an incoming call or a triggered accessory, the following happens:

- 1 The LED on the transmitter lights up in green to show that it is signalling the receiver.
- 2 The Visit LED on the receiver lights up in yellow and it starts to sound, flash or vibrate with a certain pace, called signal pattern. The transmitter and the connected accessories determine the signal pattern. The default is as follows:

| Transmitter    |               | Visit receiver |                  |           |       |
|----------------|---------------|----------------|------------------|-----------|-------|
| Source         | LED           | LED            | Sound            | Vibration | Flash |
| Landline phone | Green, top    | Yellow light   | Ring signal, low | Medium    | Yes   |
| Accessory      | Green, bottom | Green light    | Door chime, low  | Slow      | Yes   |

## Changing the signal pattern

The transmitter controls the signal pattern. Open the transmitter front cover and move the signal switches according to the tables to change it:

| Transmitter |                      | Receiver signal pattern |                       |           |       |
|-------------|----------------------|-------------------------|-----------------------|-----------|-------|
| Switch      | Source               | LED                     | Sound                 | Vibration | Flash |
|             | Landline/test button | Yellow light            | 1 x ring signal, low  | Medium    | Yes   |
|             | Mobile phone sensor  | Yellow blinks           | 2 x ring signal, high | Medium    | Yes   |
|             | Other accessory      | Green light             | 1 x door chime, low   | Slow      | Yes   |
|             | Landline/test button | Yellow light            | 1 x ring signal, low  | Medium    | Yes   |
|             | Mobile phone sensor  | Yellow blinks           | 2 x ring signal, high | Medium    | Yes   |
|             | Other accessory      | 2 x green blinks        | 2 x door chime, low   | Slow      | Yes   |
|             | Landline/test button | Yellow light            | 1 x ring signal, low  | Medium    | Yes   |
|             | Mobile phone sensor  | Yellow blinks           | 2 x ring signal, high | Medium    | Yes   |
|             | Other accessory      | 3 x yellow blinks       | 1 x ring signal, high | Medium    | Yes   |
|             | Landline/test button | Yellow light            | 1 x ring signal, low  | Medium    | Yes   |
|             | Mobile phone sensor  | Yellow blinks           | 2 x ring signal, high | Medium    | Yes   |
|             | Other accessory      | 2 x orange blinks       | Baby melody           | Fast      | Yes   |

## Transmitter

## Receiver signal pattern

| Switch                                                                            | Source               | LED               | Sound                 | Vibration   | Flash |
|-----------------------------------------------------------------------------------|----------------------|-------------------|-----------------------|-------------|-------|
|   | Landline/test button | 2 x yellow blinks | 2 x ring signal, low  | Medium ■□■□ | Yes   |
|                                                                                   | Mobile phone sensor  | Yellow light      | 1 x ring signal, low  | Medium ■□■□ | Yes   |
|                                                                                   | Other accessory      | 3 x orange blinks | Baby melody           | Fast ■■■■■■ | Yes   |
|   | Landline/test button | 2 x yellow blinks | 2 x ring signal, low  | Medium ■□■□ | Yes   |
|                                                                                   | Mobile phone sensor  | Orange blinks     | Baby melody           | Fast ■■■■■■ | Yes   |
|                                                                                   | Other accessory      | Green blinks      | 2 x door chime, high  | Slow ■□□□   | Yes   |
|   | Landline/test button | Orange blinks     | Baby melody           | Fast ■■■■■■ | Yes   |
|                                                                                   | Mobile phone sensor  | 3 x yellow blinks | 1 x ring signal, high | Medium ■□■□ | Yes   |
|                                                                                   | Other accessory      | 2 x green blinks  | 2 x door chime, low   | Slow ■□□□   | Yes   |
|  | Landline/test button | 3 x yellow blinks | 1 x ring signal, high | Medium ■□■□ | Yes   |
|                                                                                   | Mobile phone sensor  | 2 x yellow blinks | 2 x ring signal, low  | Medium ■□■□ | Yes   |
|                                                                                   | Other accessory      | Green blinks      | 2 x door chime, high  | Slow ■□□□   | Yes   |

## Transmitter

## Receiver signal pattern

| Switch                                                                               | Source               | LED               | Sound                 | Vibration   | Flash |
|--------------------------------------------------------------------------------------|----------------------|-------------------|-----------------------|-------------|-------|
|   | Landline/test button | Yellow blinks     | 2 x ring signal, high | Medium ■□■□ | Yes   |
|                                                                                      | Mobile phone sensor  | Yellow light      | 1 x ring signal, low  | Medium ■□■□ | Yes   |
|                                                                                      | Other accessory      | Green light       | 1 x door chime, low   | Slow ■□□□   | Yes   |
|   | Landline/test button | Yellow blinks     | 2 x ring signal, high | Medium ■□■□ | Yes   |
|                                                                                      | Mobile phone sensor  | Yellow light      | 1 x ring signal, low  | Medium ■□■□ | Yes   |
|                                                                                      | Other accessory      | 2 x green blinks  | 2 x door chime, low   | Slow ■□□□   | Yes   |
|   | Landline/test button | Yellow blinks     | 2 x ring signal, high | Medium ■□■□ | Yes   |
|                                                                                      | Mobile phone sensor  | Yellow light      | 1 x ring signal, low  | Medium ■□■□ | Yes   |
|                                                                                      | Other accessory      | 3 x yellow blinks | 1 x ring signal, high | Medium ■□■□ | Yes   |
|  | Landline/test button | Yellow blinks     | 2 x ring signal, high | Medium ■□■□ | Yes   |
|                                                                                      | Mobile phone sensor  | Yellow light      | 1 x ring signal, low  | Medium ■□■□ | Yes   |
|                                                                                      | Other accessory      | 2 x orange blinks | Baby melody           | Fast ■■■■■■ | Yes   |


Transmitter

Receiver signal pattern

| Switch                                                                            | Source               | LED               | Sound                 | Vibration   | Flash |
|-----------------------------------------------------------------------------------|----------------------|-------------------|-----------------------|-------------|-------|
|   | Landline/test button | 3 x yellow blinks | 1 x ring signal, high | Medium ■□■□ | Yes   |
|                                                                                   | Mobile phone sensor  | Yellow light      | 1 x ring signal, low  | Medium ■□■□ | Yes   |
|                                                                                   | Other accessory      | 3 x orange blinks | Baby melody           | Fast ■■■■■■ | Yes   |
|   | Landline/test button | Orange blinks     | Baby melody           | Fast ■■■■■■ | Yes   |
|                                                                                   | Mobile phone sensor  | 3 x yellow blinks | 1 x ring signal, high | Medium ■□■□ | Yes   |
|                                                                                   | Other accessory      | Green blinks      | 2 x door chime, high  | Slow ■□□□   | Yes   |
|   | Landline/test button | Orange/Red blink  | VMA                   | Long ■□     | Yes   |
|                                                                                   | Mobile phone sensor  | Orange/Red blink  | VMA                   | Long ■□     | Yes   |
|                                                                                   | Other accessory      | Orange/Red blink  | VMA                   | Long ■□     | Yes   |
|  | Landline/test button | Red blinks        | Fire alert            | Long ■□     | Yes   |
|                                                                                   | Mobile phone sensor  | Red blinks        | Fire alert            | Long ■□     | Yes   |
|                                                                                   | Other accessory      | Red blinks        | Fire alert            | Long ■□     | Yes   |

Trigger input

You can use both stereo and mono plugs. The 3.5 mm input is located on the top right side and the 2.5 mm input under the lid.

Mono

Trig is made by closing the two connectors

Stereo

Trig is made by:

- 1 Closing the outer and middle connector or closing the outer and inner connector
- 2 Applying 2-30VDC or 3-24VAC 5-150Hz between the middle and inner connector

Please note that when you use the 2.5 mm ext trig, the event must last for more than 2 seconds


## Changing the radio key


If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference, you need to change the radio key on all units. The radio key switches are located under the transmitter cover.

### Here is how you change the radio key:

- 1 Carefully open the front cover on the transmitter. Move any radio key switch to the up (on position) to change the radio key.
- 2 Press and hold the test/function button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 3 Press both test buttons simultaneously on the transmitter within 30 s to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Please note that all Visit units must be set to the same radio key in order to operate as a group.

## Manufacturer

Bellman & Symfon Group AB  
Södra Långebergsgatan 30  
436 32 Askim Sweden  
Phone +46 31 68 28 20  
E-mail [info@bellman.com](mailto:info@bellman.com)  
**bellman.com**

Revision: BE1431\_144MAN1.0

Date of issue: 2022-08-25

TM and © 2022 Bellman & Symfon AB.

All rights reserved.